# STANDARD HIGH SCHOOL ZZANA

# A. STANDARDISATION

1. In this experiment you are required to determine the concentration of sodium carbonate in moles per liter (molarity).

(molarity).				
You are provided with the following BA1 is a solution of sodium carbon BA2 is 0.2M hydrochloric acid so	nate			
Procedure: Pipette 20 or 25cm <sup>3</sup> of sodium car Methyl orange indicator, then titra reached.				point is
Record your results in the table be	low and repeat tit	ration until you ob	tain consist results.	<b>Results:</b>
Volume of pipette used =25.0.			cm <sup>3</sup>	
Titration number	1	2	3	
Final burette reading (cm <sup>3</sup> )	20.00	40.20	20.00	
Initial burette reading (cm <sup>3)</sup>	0.00	20.00	0.00	
Volume of BA2 used (cm <sup>3</sup> )				
Titre values used to calculate volu	me of BA2			<del></del> -
AveragevolumeofBA2				
a) Find the number of moles of hy	drochloric acid in	BA2 that reacted.		

AveragevolumeofBA2
a) Find the number of moles of hydrochloric acid in BA2 that reacted.
b) Write the equation for the reaction
write the equation for the feaction
c) Find the concentration of sodium carbonate in moles per liter (molarity)

g/ liter	•			
You are provided with the following BA1 which is a solution containing BA2 is oxalic acid solution.  Procedure; Pipette 20 or 25cm³ of BA1 into a Titrate the mixture with BA2 from Record your results in the table be	ng 2g of Sodium hy conical flask add	2drops of phenolp the end point is re	hthalein indicator ached.	ults. <b>Results:</b>
Volume of pipette used =	20.0		cm <sup>3</sup>	
Titration number	1	2	3	
Final burette reading (cm <sup>3</sup> )	11.10	21.10	31.10	
I mai barette reading (em )	11.10	21.10	31.10	
Initial burette reading (cm <sup>3)</sup>	00.00	11.10	00.00	
Volume of BA2 used (cm <sup>3</sup> )				
a).Determine the number of moles	s of sodium hydro	xide in BA1 that re	eacted	
b).Write the equation for the react	ion between BA1			
c).Determine the concentration of	oxalic acid in g po	er liter (H=1,C=12	2,O=16)	

2. In this experiment you are required to determine the concentration of oxalic acid in

# **B. PERCENTAGE PURITY AND IMPURITY**

- 1. In this experiment, you are required to determine the percentage purity and impurity in sodium carbonate.
- 1. You are provided with the following

**BA1** is made by dissolving 12.7g of the impure Sodium carbonate in one litre of solution.

**BA2** is 0.2M Hydrochloric acid solution.

#### **Procedure:**

Pipette 20 or 25cm<sup>3</sup> of BA1 into a conical flask; add 2 drops of methyl orange indicator then titrate the mixture with BA2 from the burette until the end point is reached. Record your results in the table below and repeat titration until you obtain consistent results.

Volume of pipette used =  $\dots 25.0 \dots$  cm<sup>3</sup>

Titration number	1	2	3
Final burette reading (cm <sup>3</sup> )	25.10	25.00	35.00
Initial burette reading (cm <sup>3)</sup>	00.00	00.00	10.00
Volume of BA2 used (cm <sup>3</sup> )			

Titre values used to calculate volume of BA2
AveragevolumeofBA2
a) Find the number of moles of Hydrochloric acid in BA2 that reacted
b)Write the equation for the reaction between BA1and BA2
c) Determine the concentration of sodium carbonate in moles per litre (molarity)
d) Determine the percentage purity and impurity of sodium carbonate (Na=23, C=12, O=16)

2. In this experiment, you are require carbonate.	d to determine	the percentage pur	ity and impurity ir	ı sodium
You are provided with the following; <b>BA1</b> is made by dissolving 6.4g of th <b>BA2</b> is 0.2M Hydrochloric acid solut Procedure: Pipette 20 or 25cm <sup>3</sup> of BA1 into a cor Titrate the mixture with BA2 from the below and repeat titration until you of Volume of pipette used =	ion. nical flask; add e Buirete until btain consisten	2drops of methyl of the end point is results.	orange indicator th ached. Record you	nen
Titration number	1	2	3	
Final burette reading (cm <sup>3</sup> )	20.20	40.20	20.00	
That burette reading (cm )	20.20	40.20	20.00	
Initial burette reading (cm <sup>3)</sup>	00.00	20.20	00.00	
Volume of BA2 used (cm <sup>3</sup> )				
Titre values used to calculate volume	of BA2			
AveragevolumeofBA2				
a) Find the number of moles of I	Hydrochloric a	cid inBA2 that read	cted	
	•••••			
b) Write the equation for the rea	ction between	BA1 and BA2		
			•••••	
c) Determine the concentration of	of sodium carb	onate in moles per	litre (molarity)	

d) Determine the percentage purity and impurity of sodium carbonate(Na=23,C=12,O=16)

.....

.....

# C. WATER OF CRYSTALLISATION

1. In this experiment, you are required to determine the value of n (molecules water of crystallization) in hydrated sodium carbonate.

You are provided with the following;

**BA1** is made by dissolving 28.6g of hydrated sodium carbonate Na<sub>2</sub>CO<sub>3</sub>.nH<sub>2</sub>O in one litre.

**BA2** is 0.2M Hydrochloric acid.

Procedure:

Pipette 20 or 25 cm<sup>3</sup> of BA1 into a conical flask; add 2drops of methyl orange indicator then Titrate the mixture with BA2 from the burette until the end point is reached. Record your results in the table below and repeat titration until you obtain consistent results.

Titration number	1	2	3
Final burette reading (cm <sup>3</sup> )	25.50	25.00	25.00
Initial burette reading (cm <sup>3)</sup>	00.00	00.00	00.00
Volume of BA2 used (cm <sup>3</sup> )			

Titre values used to calculate volume of BA2	
Average volume ofBA2	
a) Find the number of moles of BA2 that reacted	
b) Write the equation for the reaction between BA1 and BA2	
c) Determine the concentration of sodium carbonate in moles per litre	••••
d) Deduce the relative formula mass of Hydrated sodium carbonate Na <sub>2</sub> CO <sub>3</sub> .nH <sub>2</sub> O and Hence find the val of n in the salt (Na=23,C=12,O=16,H=1)	lue

2. In this experiment you are required to determine the value of X (molecules water of	crystallization) in
the acid	

You are provided with the following

BA1 is 0.2M sodium hydroxide solution.

BA2 is made by dissolving 12.6g of hydrated acid H<sub>2</sub>C2O<sub>4</sub>.XH<sub>2</sub>Oin one litre.

# **Procedure:**

Pipette 20 or 25cm<sup>3</sup> of BA1 into a conical flask, add 2drops of phenolphthalein indicator

Then titrate the mixture with BA2 from the burette until the end point is reached .Record your results in the table below and repeat titration until you obtain consistent results.

Titration number	1	2	3
Final burette reading (cm <sup>3</sup> )	21.00	41.00	20.00
Initial burette reading (cm <sup>3)</sup>	00.00	21.00	00.00
Volume of BA2 used (cm <sup>3</sup> )			

Titre values used to calculate volume of BA2
a) Determine the number of moles of sodium hydroxide is BA1 used
b) Write the equation for the reaction between BA1 and BA2
c) Determine the concentration of the acid in moles per dm³ (molarity)
d) Find the mass of one mole of the hydrated acid (molecular weight) and hence deduce the value of X in the acid.(H=1C=12,O=16)

# C. MOLE RATIO (stochiometry).

In this experiment you are required to determine the mole ratio between Q and P.

1. You are provided with the following;

BA1 is 0.1M solution of base Q

BA2 is 0.05M solution of an acid P

Pr	•^	^^	A.			^	
11	· U	CE	u	u	L	t	٠

**Titration number** 

Final burette reading (cm<sup>3</sup>)

Pipette 20 or 25 cm<sup>3</sup> of BA1 into a conical flask and add 2 drops of phenolphthalein indicator Then titrate the mixture with BA2 from the burette until the end point is reached. Record your results in the table below and repeat titration until you obtain consistent results.

25.00

3

35.00

1

25.50

Initial	burette reading (cm <sup>3)</sup>	00.00	00.00	10.00					
Volum	e of BA2 used (cm <sup>3</sup> )								
	alues used to calculate volume oge volume of BA2	f BA2			 				
a)	Find the number of moles of Q	inBA1 that reacte	ed						
b)	Find the number of moles of P	in BA2							
c)	determine the mole ratio between	en Q and P							
			••••••	••••••••••	••••••••••				

2. In this experiment you are required to determine the mole ratio between Q and P.
You are provided with the following;
BA1 is 0.2M solution of base R

BA2 is made by dissolving 9.8g of an acid K in one litre.

ъ	,			
Pro	red	ш	re	•

Pipette 20 or 25cm³ of BA1 into a conical flask and add 2drops of phenolphthalein indicator. Then titrate the mixture with BA2 from the burette until the end point is reached. Record your results in the table below and repeat titration until you obtain consistent results.

Volume of pipette used =  $\dots 20.0 \dots cm^3$ 

Titration number	1	2	3
Final burette reading (cm <sup>3</sup> )	20.00	40.50	20.00
Initial burette reading (cm <sup>3)</sup>	00.00	20.00	00.00
Volume of BA2 used (cm <sup>3</sup> )			

	alues used to calculate volume of BA2ge volume of BA2
a)	Find the number of moles of R in BA1 that reacted
b)	Find the number of moles of K in BA2 (RFM of $K = 98$ )
c)	Determine the mole ratio between R and K

# D. BASICITY OF AN ACID

In this experiment you are required to determine the Basicity of an acid (value of n) in H<sub>n</sub>X.

1. You are provided with the following

**BA1** is 0.1M sodium hydroxide solution

**BA2** is 0.05M solution of an acid H<sub>n</sub>X.

# **Procedure**;

Titration number	1	2	3
Final burette reading (cm <sup>3</sup> )	25.70	25.00	25.00
Initial burette reading (cm <sup>3)</sup>	00.00	00.00	00.00
Volume of BA2 used (cm <sup>3</sup> )			

	alues used to calculate volume of BA2ge volume of BA2
a)	Find the number of moles of Sodium hydroxide in BA1 that reacted
b)	Find the number of moles of the acid that reacted.
c)	Find the value of <b>n</b> in the acid.

BA1 is 0.2M BA2 is 0.1M	ided with the following sodium hydroxide so solution of an acid H	olution					
	25cm <sup>3</sup> of BA1 into a						
	n titrate it with BA2 f e table below and repo		-	-	your		
Volume of pi	pette used =	20.0		cm <sup>3</sup>			
Titration nu	mber	1	2	3			
Final burette	reading (cm <sup>3</sup> )	42.00	40.00	40.00			
Initial burette	e reading (cm <sup>3)</sup>	00.00	00.00	00.00			
Volume of B.	A2 used (cm <sup>3</sup> )						
Titre values v	used to calculate volume of BA2	me of BA2					
a) Find t	the number of moles	of Sodium hydrox	ide in BA1 that re	acted			
	b) Find the number of moles of the acid that reacted.						
	the value of n in the a						
•••••							

2. In this experiment you are required to determine the Basicity of an acid (value of n) in  $H_nX$ .

# E. DETERMINATION OF RELATIVE ATOMIC MASS

1. In this experiment you are required to determine the relative atomic mass of X in one mole of l	<b>1.</b> ]	In th	nis ex	periment	you are	required	to d	letermine	the	relative	atomic	mass	of X ir	one	mole	of F	$\zeta_2$	X
--	-------------	-------	--------	----------	---------	----------	------	-----------	-----	----------	--------	------	---------	-----	------	------	-----------	---

You are provided with the following:

**BA1** which is a solution containing 1.725g of  $K_2X$  in  $250cm^3$  of solution.

**BA2** which is 0.1M sulphuric acid.

Sulphuric acid reacts with K<sub>2</sub>X in the ratio 1:1

# **Procedure:**

Pipette 25 cm<sup>3</sup> or 20 cm<sup>3</sup> of BA1 into a conical flask, add 2 or 3 drops of methyl orange indicator and titrate with BA2 from the burette. Repeat the titration until you obtain consistent results. Record you obtain consistent results. Record your results in the table below.

Final burette reading/cm <sup>3</sup>	10.00	21.00	31.00
Initial burette reading/cm <sup>3</sup>	0.00	10.00	21.00
Volume of BA2 used/cm <sup>3</sup>			

Volumes of BA2 used in the calculating the average volume
Average volume of BA2
a)Calculate number of moles of BA2 that reacted.
b) Determine the;
i) molarity of BA1
ii) the relative formula mass of $K_2X$ and hence the mass of X in one mole of $K_2X$ (K=39)

2 - \$7.		Съд		шсо					
2. <b>You are required to deter</b> You are provided with the follows:		mass of M in the	metal carbonate M	IHCU3					
<b>BA1; which</b> is 0.2M hydrochlo	· ·								
BA2; which is a solution conta		3 of a motal hydro	gan aarbanata MUC	Λ.					
Procedure:	ining 6.4 g/300cm	of a metal flydio	igen carbonate winc	$O_3$					
	:	1 11 2 1 f							
Pipette 25.0 cm <sup>3</sup> or 20 cm <sup>3</sup> of BA2		-	•						
the burette. Repeat the titration unt	•		•	e table below.					
Volume of pipette used =	.25.0		cm <sup>3</sup>						
Titration number	1	2	3						
Final burette reading (cm <sup>3</sup> )	25.00	27.00	25.20						
-									
Initial burette reading (cm <sup>3)</sup>	0.00	0.00	0.00						
W. L. (CDA) 1 ( 3)									
Volume of BA2 used (cm <sup>3</sup> )									
Titre values used to calculate volum	ne of BA2								
Average volume of BA2									
Write the equation for the react	ion between hydro	ochloric acid and N	$MHCO_3$						
a) Calculate the number of mo	oles of BA1 that re	acted.							
	•••••	•••••							
	•••••	•••••	•••••						
1) 5	•••••	••••••							
b) Determine the;									
i) molarity of BA2									
	•••••	•••••	•••••						
	•••••	•••••	•••••						
		•••••	•••••						
	•••••	•••••	•••••	•••••					

ii) the relative formula mass of the hydrogen carbonate and hence the atomic mass of M (H=1,C=12,O=16)

.....